CAWO DIGITIZER CR2000ME

High Resolution CR-System

CAWO DIGITIZER CR2000ME

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... is a modern, high resolution scanner for imaging plates. It transfers x-ray information saved on the imaging plate into a digital signal. The CONAXX Software (included in the standard delivery) operates the CR-Scanner and routes the images automatically to a PACS. In case you do not have a PACS yet, our PROMIS medical imaging software would be a favorable alternative.

Rapid imaging for superior quality

Due to the fast scanning process, x-ray images are provided quickly for a reliable diagnosis. Misdosage and as a result overexposure or underexposure are almost impossible with this CR-technology due to the high exposure tolerance.


Small investment

The worthwhile investment brings the immediate benefit of no longer needing to work with chemicals, the disposal of chemistry and the cleaning of a processor. Expensive printouts on film are not required, diagnosis can be done fast and reliably on a diagnostic monitor. The imaging plates can be re-used and you can still continue to use your existing x-ray unit.

Easy changeover for the radiographer because the handling is comparable to that of films.

Reliable diagnosis on the monitor (even with remote access). Printouts are not needed.

Highest precision due to state of the art technology and certified quality management.

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Human medicine

In human medicine the CAWO CR2000ME CR-System could be used either as independent modality, integrated in a PACS environment via CONAXX or as a stand alone solution for a private practitioner. Through the operation of the CR-Scanner and image distribution with CONAXX Software the CR2000ME can be used very flexibly. The images are just routed to an existing PACS. If there is no PACS available yet, a system for diagnosis and archiving of images can be customized according to the needs of the practice or clinic with our modular structured PROMIS medical imaging software. Connections to patient/practice administration software can easily be realized, too.

Veterinary medicine

Veterinarians appreciate the advantages of the CAWO CR2000ME CR-System as well. Especially when it comes to mobile x-rays. The uncomplicated handling and the high mobility (transportation case is available as an accessory) simplify the workflow when x-raying large animals. A diagnosis can be made on site immediately after taking the x-ray (e.g. on a laptop monitor) and further treatment can be initiated right away.

Best image quality for all standard x-ray applications in human medicine.

CR2000ME offers affordable CR-Technology for the ambitious veterinarian.

Due to its compact design the CR2000ME is suitable for **mobile** applications.

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CR2000ME – revolutionary PCS-technology

The revolutionary PCS-technology (Photon-Collecting-System) ensures extraordinary low image noise by scanning the imaging plate with almost 100% efficiency of active radiation. Due to the very high resolution of more than 5 line pairs per millimetre an exact diagnosis is ensured. An outstanding feature is the greyscale resolution of 16-Bit. The CR2000ME CR-System uses a centric laser beam, which is guided towards the image plate by means of a rotating pentagonal prism. The high rotation speed of this prism in connection with the centrally located digital sensor (PMT) that collects the luminescence radiation (emitted by the IP after being activated with the laser) and the linear drive that gently feeds the IP through the scanner ensure the outstanding resolution. Immediately after scanning the image data is automatically erased from the imaging plate. The IP can be used for the next x-ray right away, erasing as an additional step is not required.

Scanning in a matter of seconds with patented laser technology and immediate erasing of image data from the IP.

Digital Imaging Plates CAWO DIGIT

CAWO is one of the most important manufacturers of digital imaging plates. We offer a wide range of products for many markets and applications.

CAWO is production partner for a number of important global players on the digital imaging field.

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Simple application. That's the way it works!

1. Exposure of the CAWO DIGIT imaging plate.
2. Parameter selection and CR2000ME operation with CONAXX.
3. Remove the imaging plate (in protection sleeve) from the cassette.
4. Placement of the imaging plate, the scanning process starts automatically.
5. Diagnosis on the screen, e.g. with PROMIS PACS.

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

Simple application, simple workflow with CAWO CR2000ME and CONAXX

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The virtual waiting room provides an overview of all upcoming x-ray jobs.

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Body part selection supports creation of individual scan jobs.

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The patient's job management screen provides a status-overview of the individual jobs.

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In the General Job Management jobs can be selected according to their status.

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In the preview the image can be processed and optimized.

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At the end of the workflow further steps are defined (e.g. routing to PROMIS).

CONAXX Acquisition System

The CONAXX Acquisition System is included in the standard delivery and operates the complete CR2000ME CR-System. On a touch-screen (optional accessory) the software guides you in a user-friendly way through the complete workflow according to the process in a practice or clinic. The CONAXX Acquisition System consists of the Image Acquisition Module as well as DICOM Store, GDT Data Input/Response and E-mail Modules. With CONAXX multi-clinic applications, e.g. when doctors share modalities, can easily be realized (optional Multi-Clinic-Module). In addition it offers the following functions:

Patient acceptance / identification

- manual entering of new patients
- acceptance from patient/practice administration systems via GDT
- automatic work list acceptance (SCU) including transfer of jobs (DICOM Input Module, optional)
- DICOM Query (part of optional DICOM Input Module)

Creation of x-ray jobs

- creation of individual scan-jobs with body part selection, incl. allocation of correct scan mode
- general and patient orientated job management

Image preview and optimization

- real time image preview
- “crop”-function (cutting of images for data reduction)
- comments, turn, flip, mirror, L/R-signs and annotations
- automatic window levelling and visual histogram
- brightness and contrast
- selection of several standard filters
- x-ray journal with automatic transfer of dose-area-product (optional)

Data transfer and archiving

- creation of DICOM-file with header information
- feedback to worklist-server (DICOM Input Module) or patient/practice administration system (GDT Data Input/Response Module incl.)
- image export as DICOM or in other formats (JPG, TIF, BMP)
- DICOM Store - transfer of images to a DICOM compliant PACS (e.g. PROMIS)
- internal “safety-archive” with automatic memory adjustment
- printing images on Windows printer or sending images by e-mail

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Tel : 086 338 1481,662 363 8043

Technical Data:

Resolution:	varies depending on size and type of imaging plate
Smallest physical pixel size:	12,5 µm (respectively 40 Lp/mm)
Max. infeed width:	35 cm (14")
Dimensions (H/W/D):	39 cm x 38 cm x 52 cm (15 1/8" x 15" x 20 1/2")
Weight:	20,4 kg
Installation:	table top
Noise:	approx. 49 dB(A)
Temperature range:	+10 up to +30 °C (+50 up to +86 °F)
Power supply:	100-240 V / 50-60 Hz
Electrical connected load:	120 W
Laser protection class:	I (EN60825-1)
Further detail:	integrated printer port

Interfaces:

Integration in existing PACS / RIS / HIS
Takeover of worklists (DICOM)
Connection to existing system via GDT/BDT

Standard delivery includes:

CONAXX software, including the following modules:
Standard, DICOM Store, GDT Data Input/Response & E-mail
2 imaging plates 18 x 24 cm / 24 x 30 cm, at buyer's option, (8" x 10" / 10" x 12")
2 imaging plates 35 x 35 cm / 35 x 43 cm, at buyer's option, (14" x 14" / 14" x 17")
4 protection sleeves (according to chosen formats)
Connecting cable (3m) for USB and power supply

Additional requirements

PC
Diagnostic monitor
Imaging storage
Film cassettes without intensifying screens

Minimal PC-system requirements:

1,4 GHz processor
512 MB main memory (RAM)
20 GB free hard disk space (recommendation >160 GB)
Operating system Windows 2000™ or Windows XP™
USB 2.0
optional: TCP/IP-connection to network

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Industrial Instrument Co.,Ltd.